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# Product Specification for Cartridge System for Pen-injectors

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### 文件批准 Document Approval

	部门 Department	姓名 Name	签字 Signature	日期 Date
起草人 Prepared by	研发部 Research and Development Department	赵自强		
审核人 Reviewed by	研发部 Research and Development Department	张鑫		
批准人 Approved by	质量部 Quality Department	孟壮		

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	部门 Department	签字 Signature	日期 Date
审核人 Reviewed by			
审核人 Reviewed by			
批准人 Approved by			

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Production address: No. 1000, Qihuang 1st Road, Shuangliu District, Chengdu, Sichuan Province

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This document applies to 3ml Cartridge system for pen-injectors (referred to as cartridge) manufactured by Sichuan Tofflon Pharmaceutical Packaging Co., Ltd.

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Cartridge system for pen-injectors are single-use, Ready to Use packaging materials for injections.

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The filling and closure of the cartridge is not completed by Tofflon and is not Tofflon's responsibility

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





To protect the product from contamination and damage, this product is packaged in nests, tubs and double steribags, and the process of capping and primary packaging are carried out in Class A clean environment in Class C background. Tub, Tyvek Lid, and steribag form the sterile barrier system. After the primary packaging, the product is loaded into a five-layer corrugated box, and after the product pallet is stacked, all corners are protected by stretch film, so as to better ensure that the product is not damaged and polluted during transportation and storage.

本产品标识信息如下：Packaging identification may include the following information:

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Note: "+" indicates that the information is identified, and "-" means that the information is not identified.

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This product should be stored at room temperature (10°C-30°C).

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Tofflon conducts long-term tests and accelerated tests on this product in accordance with T/CNPPA 3017-2021 "Guidelines for Self-Stability Research of Plastic and Rubber Pharmaceutical Packaging Materials" to evaluate the self-stability of the product. The shelf life of this product is tentatively set at 3 years.

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Tofflon sets product quality standards and inspects and releases the product in accordance with the general principles of the ChP, YBB00132004-2015, ISO 13926-1, USP<660>, EP 3.2.1 and ISO 21881.

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The compatibility study of this product is conducted in accordance with the "Technical Guidelines for the Compatibility Study of Chemical Drug Injections and Medicinal Glass Packaging Containers (Trial)", "Technical Guidelines for the Compatibility Study of Chemical Drugs and Elastomer Lids (Trial)", and ICH Q3D.

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Tofflon establishes, maintains and continuously improves our quality management system in accordance with regulations and standards such as "Good Manufacturing Practice for Pharmaceutical Manufacturing", "Good Manufacturing Practice for Pharmaceutical Packaging Materials" and ISO 15378 *Primary packaging materials for medicinal products-Particular requirements for the application of ISO 9001:2015, with reference to good manufacturing practice(GMP)* to ensure that we continue to provide products and services that meet the requirements and have controllable quality.

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This product complies with the following regulations and standards:

- 中国药典《笔式注射器用卡式瓶系统指导原则》；
- ChP Guideline on Cartridge system for pen-injectors
- 中国药典 9622《药用玻璃容器指导原则》；
- ChP 9622 Guideline on Glass Containers for Pharmaceutical Packaging
- 中国药典《药品包装用橡胶密封件指导原则》
- ChP 9623 Guideline on Rubber Closures for Pharmaceutical Packaging
- YBB00132004-2015《笔式注射器用硼硅玻璃套筒》；
- YBB00132004-2015 Borosilicate Glass Barrels for Pen-injector
- YBB00162004-2015《笔式注射器用溴化丁基橡胶活塞和垫片》
- YBB00162004-2015 Bromobutyl Rubber Plunger and Liner for Pen-injector
- ISO 13926 Pen systems
- USP<660> Containers-Glass
- EP 3.2.1 Glass containers for pharmaceutical use

- ISO 21881 Sterile packaged ready for filling glass cartridges

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Tofflon takes the following measures to control the quality of this product from the design:

- 采购已在上市药品中使用多年的原材料；
- Procurement of raw materials that have been used in marketed drugs for many years
- 套筒经清洗后，通过隧道烘箱进行烘干并去除热原；
- After the barrel is washed, it is dried in a tunnel oven and pyrogen is removed
- 产品清洗烘干后的轧盖至内包装工序在 C 级背景下的 A 级洁净环境中进行，最大限度地控制产品微粒、细菌内毒素和微生物污染；
- After drying process, the capping and primary packaging is carried out in a Class A clean environment in Class C background, so as to control the product particles, bacterial endotoxins and microbial contamination to the greatest extent
- 生产线上配置了多个在线外观检测装置，第一时间剔除不良品；
- The production line is equipped with multiple appearance inspection devices to reject defective products at the first time;
- 采用全自动运线，杜绝人为差错。
- Automatic transportation line is adopted to eliminate human error.

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Tofflon formulate the inspection rules and Acceptance Quality Limit of the product inspection items in accordance with the ChP "Guiding Principles for the Inspection Rules of Pharmaceutical Packaging Materials" and GB/T 2828 (idt.: ISO 2859), as shown in the following table.

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## 活塞要求 Requirements of Plunger

检验项目 Items	接受标准 Acceptance criteria	缺陷等级 分类 Defect grade classification	是否批 检验 Batch inspection?	检验水 平 Inspection level	接收质量 限 (AQL)
外观 Appearance	本品表面色泽应均匀，不得有污点、杂质、气泡、裂纹、缺胶、粗糙、胶丝、胶屑、海绵、毛边；不得有除边造成的缺陷或锯齿现象；不得有模具造成的明显痕迹 The color of this product should be uniform, and there should be no stains, impurities, bubbles, cracks, lack of glue, roughness, rubber wire, rubber crumbs, sponges, and burrs; There shall be no defects or jagged phenomena caused by edge removal; There must be no obvious marks during molding.	随机性缺陷 B 类 Random Defect B	是 Yes	I	0.65
尺寸 Dimension	应符合质量标准规定 Should meet the requirements of quality standard	随机性缺陷 A 类 Random Defect A	是 Yes	S-2	0.4
红外光谱 Infrared spectrum	红外图谱均应与典型图谱基本一致 The infrared spectrum should be basically consistent with the typical spectrum	系统性缺陷 Systemic Defect	否 No	N/A	N/A
灰分 Ash	标称值±2% Nominal value±2%	系统性缺陷 Systemic Defect	否 No	N/A	N/A
挥发性硫化物 Volatile sulfides	不得深于标准硫斑 Not be deeper than standard sulfur spots (1.0µg/cm <sup>2</sup> )	系统性缺陷 Systemic Defect	否 No	N/A	N/A
不溶性微粒 Insoluble particles	≥10µm 的微粒数: ≤60 粒/ml; ≥25µm 的微粒数: ≤6 粒/ml ≥10µm particle NMT 60 /ml ≥25µm particle NMT 6 /ml	随机性缺陷 A 类 Random Defect A	是 Yes	N/A	N/A
细菌内毒素 BET	<0.25EU/ml	随机性缺陷 A 类 Random Defect A	是 Yes	N/A	N/A
无菌 Sterility	应符合规定 Be compliance	随机性缺陷 A 类 Random Defect A	是 Yes	N/A	N/A

## 活塞溶出物要求 Extractables test of Plunger

检验项目 Items	接受标准 Acceptance criteria	缺陷等级 分类 Defect grade classification	是否批 检验 Batch inspection?	检验水 平 Inspection level	接收质量 限 (AQL)
澄 清 度 与 颜 色 Clarity and color	供试液应澄清无色。如显浑浊，与 0.5 号浊度标准液比较，不得更浓。如显色，与黄绿色 5 号标准比色液比较，不得更深。 The test solution should be clear, any turbidity produced should be not more intense than turbidity standard solution No.0.5.The test solution should be colorless, any color observed should be not more intense than yellow-green No. 5 standard colorimetric solution.	N/A	N/A	N/A	N/A
pH 变化值 pH change value	供试液与空白液 pH 值之差不得过 1.0 The difference between the pH value of the test solution and the blank solution shall not exceed 1.0				
吸光度 Absorbance	在波长 220~360nm 范围内最大吸光度应≤0.1 The maximum absorption value between 220 and 360nm shall not exceed 0.1.				
易氧化物 Oxidizable substance	供试液与对照液消耗硫代硫酸钠滴定液 (0.01mol/L) 体积之差不得过 3.0mL The difference in volume between the consumed sodium thiosulfate titrant (0.01mol/L) of the test solution and the control solution should not exceed 3.0mL.				
不挥发物 Nonvolatile matter	供试液与空白液遗留残渣之差不得过 2.0mg The difference between the residue left by the test solution and the blank solution should not exceed 2.0mg				
电导率 Electrical conductivity	≤20.0μS/cm				
铵离子 Ammonium ion	与对照液目视比色，不得更深 (0.0002%) Compared with the control solution visually, it should not be darker (0.0002%)				
重金属 Heavy metals	供试液中重金属含量不得过百万分之一 ≤1ppm				
锌离子 Zinc ion	供试液应不得深于对照液。(0.0003%) Compared with the control solution				

检验项目 Items	接受标准 Acceptance criteria	缺陷等级 分类 Defect grade classification	是否批 检验 Batch inspection?	检验水 平 Inspection level	接收质量 限 (AQL)
	visually, it should not be darker (0.0003%)				

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



**Tofflon** 东富龙

此图纸版权属于东富龙所有, 对未经许可擅自拷贝、使用者, 本公司保留追究其法律责任的权利。  
THIS DRAWING IS COPYRIGHT. IT MUST NOT BE COPIED, REPRODUCED OR OTHERWISE MADE USE OF  
WITHOUT THE TOFFLON'S WRITTEN PERMISSION.

标题/TITLE:

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页/SHEET: 1 共/OF: 1 图幅/SIZE: A4

图号/DWG: TZ-TFL-1044-01

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01

客户名/CUSTOMER:

项目号/PROJECT NUMBER:

设计/DESIGNED: wanglei

日期/DATE: 2023/06/12

质量/MASS:

未注明公差表  
TOLERANCED UNLESS OTHERWISE NOTED

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

**Tofflon** 东富龙

此图纸版权属于东富龙所有，对未经许可擅自拷贝、使用者，本公司保留追究其法律责任的权利。  
THIS DRAWING IS COPYRIGHT. IT MUST NOT BE COPIED, REPRODUCED OR OTHERWISE MADE USE OF  
WITHOUT THE TOFFLON'S WRITTEN PERMISSION.

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



**This product is composed of two parts: rubber gasket and aluminum cover, and the gasket is a double-layer gasket, which is polyisoprene rubber and halobutyl rubber.**

2. 内外表面应清洁，无残留溶剂、毛刺和损伤 **The inner and outer surfaces should be clean and free of residual solvents, burrs, and damage.**

<b>Tofflon</b> 东富龙		此图纸版权属于东富龙所有，对未经许可擅自拷贝、使用者，本公司保留追究其法律责任的权利。 THIS DRAWING IS COPYRIGHT. IT MUST NOT BE COPIED, REPRODUCED OR OTHERWISE MADE USE OF WITHOUT THE TOFFLON'S WRITTEN PERMISSION.							
标题/TITLE: <b>笔式注射器用含垫片铝盖</b> Aluminum Cap with Gasket for Pen-injector		页/SHEET: 1 共/OF: 1	图幅/SIZE: A4_CN						
客户名/CUSTOMER:		项目号/PROJECT NUMBER:	图号/DWG: TZ-TFL-2035-01						
设计/DESIGNED: ZX	日期/DATE: 2025.01.10	质量/MASS:	版本/REV: 01						
校核/CHECKED: WJY	日期/DATE: 2025.01.10	比例/SCALE: <b>8:1</b>	未注明公差表 TOLERANCED UNLESS OTHERWISE NOTED						
审核/APPROVED: MZ	日期/DATE: 2025.01.10	材料/MATERIAL: <b>Al+Rubber</b> 铝/橡胶	<table border="1"> <tr> <td>MACHINING 机加工</td> <td>WELDING 焊接</td> <td>SHEET METAL 钣金</td> </tr> <tr> <td>GB/T 1804-m GB/T 1184-K</td> <td>GB/T 19804-B GB/T 19804-F</td> <td>GB/T 15055-m</td> </tr> </table>	MACHINING 机加工	WELDING 焊接	SHEET METAL 钣金	GB/T 1804-m GB/T 1184-K	GB/T 19804-B GB/T 19804-F	GB/T 15055-m
MACHINING 机加工	WELDING 焊接	SHEET METAL 钣金							
GB/T 1804-m GB/T 1184-K	GB/T 19804-B GB/T 19804-F	GB/T 15055-m							

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The dimensions of this product comply with ISO 13926-2:2017

<b>Tofflon</b> 东富龙		此图纸版权属于东富龙所有，对未经许可擅自拷贝、使用者，本公司保留追究其法律责任的权利。 THIS DRAWING IS COPYRIGHT. IT MUST NOT BE COPIED, REPRODUCED OR OTHERWISE MADE USE OF WITHOUT THE TOFFLON'S WRITTEN PERMISSION.		
标题/TITLE: BROMOBUTYL RUBBER PLUNGERS FOR PEN-INJECTORS(3ML)		页/SHEET: 1 共/OF: 1	图幅/SIZE: A4_CN	
客户名/CUSTOMER:		图号/DWG: TZ-TFL-1048-01		版本/REV: 01
设计/DESIGNED: ZX		日期/DATE: 2024.04.12	项目号/PROJECT NUMBER:	
校核/CHECKED:		日期/DATE:	质量/MASS:	
审核/APPROVED: SYT		日期/DATE: 2024.04.12	比例/SCALE: <b>5:1</b>	
材料/MATERIAL: BROMOBUTYL RUBBER		未注明公差表 TOLERANCED UNLESS OTHERWISE NOTED		
		MACHINING 机加工	WELDING 焊接	SHEET METAL 钣金
		GB/T 1804-m GB/T 1184-K	GB/T 19804-B GB/T 19804-F	GB/T 15055-m